SPECIFICATION

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RF-MTQAB7N03-B1

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Mass Product 量产供货

Contents 目录

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1. Description 产品介绍

1.1 General Description 产品描述

Notes 备注: The product consists of LED module (red light, blue light, IR), mask and control panel
产品由 LED 模组 (红光、蓝光、IR) 和面罩以及控制板组成的 LED 面部美容仪

1.2 Features 产品特征

- ▶ Appearance of private module design, fashion and beautiful;
外形私模组设计，时尚美观;
- ▶ According to the facial uniform design, so that every inch of skin feel light;
根据面部均匀设计，使每一寸肌肤感受到光;
- ▶ Low voltage DC3.7V constant voltage driving to ensure safe and reliable running ;
低压 3.7V 工作，确保运行安全可靠；
- ▶ Visual design, beauty while working;
可视设计，边工作边美容;
- ▶ Built-in battery, easy to carry ;
内置电池，方便携带.
- ▶ Touch switch, switch at will.
触控开关，随心切换.

1.3 Application 产品应用

- ▶ Facial beauty app

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1.4 Features Introduction 产品功能

▶ Red light, 650-660nm, promotes blood circulation, accelerates metabolism, promotes collagen production under the skin, and reduces fine lines and wrinkles;

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▶ Blue light 405-420nm, improve the pore environment, improve acne, reduce the frequency of acne;

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
▶ Powder light +IR (red light & Blue light & IR) can improve the alternating function of cell oxygen, improve microenvironment, fade pigment, soften wrinkles and moisturize skin ;

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
Fig.3-2 Title 图名

5. Handling Precautions 产品使用注意事项

5.1 Product storage and maintenance 产品存储和维护

▶ Do not expose the product to rain, humidity and high temperature.

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▶ Open flame sources, such as lit candles, should not be placed on the device

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▶ Careful supervision is required to prevent children from using this product.

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▶ If not used for a long time, charge the battery at least once every 3 weeks to maintain the battery life

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▶ Do not poke this product with sharp metal objects such as keys, hairpins, necklaces or nails.

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▶ Do not put this product in the oven, dishwasher, washing machine

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5.2 Handling Precautions 产品使用注意事项

Products for facial care, not recommended for body use.

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Do not use this device for skin allergies, wounds, severe acne, burns, and bacterial skin infections.

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Please consult a dermatologist before using this device if you are being treated for a skin disease or have a skin condition.

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Not recommended for pregnant, menstruating, lactating women (unstable hormonal balance may cause temporary skin problems).

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Do not use for people who are unable to express themselves clearly (including babies).

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People with light sensitivity should not use this product if they are allergic to ultraviolet light.

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Do not use this product in the rain, damp (bathtub, shower, sink), high temperature (open flame source, oven) environment.

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Do not disassemble without permission, there is no user repairable parts inside, related maintenance, please consult qualified maintenance personnel.

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REFOND

Declare 申明

This specification is written both in English and in Chinese and the latter is formal.
产品规格书以中英文方式书写，若有冲突以中文版本为准。